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1. Instructions for Use (EN) | Explorer System | Socket

Valid from: 01.02.2024

- Read this document carefully before using the product.
- Observe the safety instructions to avoid injuries and product damage.
- Instruct the user in the proper and safe use of the product.
- Contact the manufacturer if you have any questions about the product (e.g., commissioning, use, maintenance, unexpected operation or incidents). You will find the contact details on the back.
- Keep this document in a safe place.

These instructions for use provide you with important information on the use, adjustment and handling of the Explorer Prosthesis System.

40. Product description

System

The Explorer Prosthesis System is a modular passive functional forearm prosthesis, hereinafter referred to as the Explorer System. The Explorer System consists of a socket design and several fitting parts, hereinafter referred to as Explorer Modules. Various modules are available for the Explorer System. The modules support the user in various activities as an aid.

Socket

The Explorer socket is custom-made by a specialist for each user and is therefore a special design. It is used with a socket ring and passive functional hands (fitting parts) that serve as accessories of a medical device.

Wrist "Coupling"

The coupling is not an actual part, but refers to the interface through which the socket is connected to one of the modules at a time. The socket is equipped with a nut, optionally with a metric (M12x1.5) or inch (1/2"-20) thread, and a socket ring. The modules are equipped with a screw, optionally with a metric (M12x1.5) or inch (1/2"-20) thread. The socket ring is used to position the module used more easily in the desired orientation.

Possible combinations

The system has been tested for possible combinations of the Explorer socket with the various Explorer modules. See chapter 12 for more information on the load capacity forces of the modules.


41. Use

Intended Use

The macu4 Explorer socket is a custom-made prosthesis socket intended for prosthetic replacement of the forearm.

Indication

The Explorer System can be used for unilateral amputees from amputation height transradially or with dysmelia in forearm fittings (according to ISO 8548-2:2020) as well as for categories 2 to 5 as shown in the diagram below.

Graph 1. Categories of forearm deficits

The system is suitable for persons 2 years and older. An individual evaluation by the specialist is necessary to assess whether the user is physically and mentally in the condition to use the system properly.

Contraindication

- The Explorer System is not indicated for individuals with bilateral forearm deficits.
- The Explorer System is not indicated for individuals who have musculoskeletal limitations (e.g., stiffness in the shoulder joint).
- The Explorer System is not indicated for individuals with fresh wounds, allergies or sensitive skin on the arm.
- The current design of the Explorer socket is not suitable for arm deformities that include, for example, a wrist and one or more residual fingers.

Operating conditions

The Explorer System was developed especially for everyday, sports and leisure activities with low to moderate stress. The system must not be used for extraordinary activities. These extraordinary activities include, for example, sports with excessive coupling and/or impact loads (push-ups, downhill, mountain biking, etc.) or extreme sports (free climbing, paragliding, etc.). In addition, the system should not be used for driving motor vehicles, driving heavy equipment (e.g., construction machinery), operating industrial machinery or operating motor-driven work equipment.

Operating conditions -- socket

The socket is pulled over a prosthetic liner, for example, then tightened using the lacing system, and only then is the hand coupled to the socket using the healthy arm, for example.

42. Warnings

Notes on unpacking the supplied components

Only open the packaging at the intended points with care so that the parts inside are not damaged. Inspect the delivered components and make sure that they are undamaged. If the components are damaged, please follow the steps in chapter 11.

Notes on commissioning the socket

Please familiarize yourself with the use of the socket before using it for the first time. Make sure that you always use the socket ring together with the modules when you use the modules together with the Explorer socket. The socket ring allows easier fixing of a module in the desired orientation.

Instructions for putting on the socket

If your arm has changed in shape, circumference and/or length since custom fabrication, the fit of the socket may no longer be optimal. Make sure that the socket encloses your arm appropriately. See chapter 6 for more information.

The socket must be positioned in a specific orientation on the arm before you use it. Make sure that the orientation is correct. Refer to chapter 5.

Instructions for checking before use

Inspect the modules, socket and other accessories before each use and make sure that they do not show any visible damage. Improper handling of the socket can restrict its function or render it inoperable. Please also observe the information on the mechanical load capacity of the socket (see chapter 11). The 3D printed component of the socket is made of polyamide 12. The mechanical stability of the socket may be impaired if it is used and/or stored outside the permissible temperature range. If possible, UV radiation should be avoided during use, but absolutely prevented during storage so as not to accelerate the aging process.

Notes on the use of the socket

Use the socket with a protective insert, prosthetic sock or liner. Do not use the socket if you feel that it no longer fits comfortably in terms of length, circumference or due to pressure points.

Notes on the correct orientation of the socket

Make sure that the button of the lacing system always points laterally (outside of your arm). See also chapter 6.

Notes on the lacing system

During an activity, be careful not to move the button of the lacing system too close to objects that could snag the button, for example, a net or a more tightly woven rope construction. Do not use the socket if the lacing system is not functional because the laces cannot be tightened or loosened.

Before putting on the upper for the first time, check that the lacing system is working properly. This also includes the complete tightening of the lacing. Even from this state, the system should now be able to be loosened again by first turning the head in the opposite direction until a click is heard and then pulling the two protruding ropes one after the other.

Notes on the secure fit of the socket

Make sure that both the lacing system and the locking mechanism on the back of the upper arm are used or closed by you. Both serve to ensure that the socket sits firmly and securely on your arm.

Notes on the condylar braces

During an activity, be careful not to move the condylar braces too close to objects that they could catch on. For example, use caution if you move the arm on which you wear the brace

Explorer System – Socket

very close to a net or a more tightly woven rope construction back and forth (sideways or forward/backward).

Hints for taking off the socket

Make sure that both the lacing system and the locking mechanism on the back of the upper arm are opened or loosened by you. Both serve to ensure that you can easily remove the upper.

Notes on the care of the socket

Make sure that you always clean the socket after use. Remove visible dirt particles so that the socket is not restricted in its function in the area of the condylar clips, the lacing system and the coupling. See chapter 6 for more information.

Notes on ambient conditions

Make sure that you store or use the socket in accordance with the information on temperature, UV radiation and humidity given in chapters 10 and 11.


Notes on magnets in product

The socket ring contains mini magnets (5x5x5mm). Make sure that you keep the socket ring at least 5cm away from e.g. pacemakers or defibrillators.

43. Application of the socket

Before using the Explorer socket for the first time, familiarize yourself with the instructions for use. The following chapters explain how to use the Explorer socket.

Socket structure

Velcro closure

Condylar braces

Proximal end of the socket

Button of the lacing system

Distal end of the socket


Socket ring

Prepare lacing system


Check that the lacing is loose enough.

To loosen the lacing, turn the knob counterclockwise about 10 degrees until you hear a soft click.


Now the cords can be loosened by gently pulling them with your fingers.


Note that the socket is open enough so that you can easily insert the arm. Press the sides of the socket slightly outwards. Increase the width as needed.

Open Velcro



Take your macu4 socket and place it between your knees, for example.

Open the Velcro fastener of the condylar braces and loosen the Velcro fastener



Tighten the socket

Put on your prosthetic sock or inlet. Carefully insert the arm into the socket. The center of the laces should be on the top of the arm. The lace button should sit on the outside of your arm (not facing the center of your body)

Insert the arm into the socket until the clips encircle the epicondyles of the elbow on both sides and fit securely. The end of the condylar clips should face the back of your arm (the Velcro runs along the back of the arm).


Lacing system and Velcro closure

When the arm is comfortably positioned in the socket, turn the knob clockwise to tighten the lacing system. Turn until the socket fits snugly on the arm. The fit of the socket on the arm should be snug but not uncomfortable.


Avoid that a too tight lacing squeezes the tissue of your arm in the area of the crook of the arm. Use the Velcro to close the socket at the back and above your elbow. This will give a tighter fit.

Final control


- Check the fit of the socket before you start any activity.
- Make sure it is well positioned and secured by trying to move the socket gently.
- Note that the socket should remain firmly in position.
- If the upper is not tight and can be moved, recheck the lacing system and Velcro closure.


Take off the socket


Open the Velcro fastener.

To loosen the lacing (e.g., Atop System), turn the knob approximately 10 degrees counterclockwise until you hear a soft click.

Now the cords can be loosened by gently pulling them with your fingers.

Note that the socket is open enough so that you can easily pull the arm out. If necessary, press the sides of the socket slightly outwards.

44. Care of the socket

Dirt particles in the spaces between the products should be removed with a soft-bristle brush. This pre-cleaning of the products should be done without the use of water. Product surfaces should be cleaned with a damp cloth (non-scratching). Cold or lukewarm water should be used for this purpose.

45. Adverse event reporting

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the member state where the user and/or patient is established. Please use the form Please use the Product Experience Report (PER)

https://25462115.fs1.hubspotusercontent-eu1.net/hubfs/25462115/mac4_AG_Explorer/Product_Experience_Report_EN.pdf and send it to by email to support@macu4.com.

46. Legal notice

Limitation or exclusion of liability

macu4 AG, hereinafter referred to as Macu, is only liable if the specified handling and processing instructions as well as the care instructions and maintenance intervals for the

product are observed. Macu expressly points out that this product may only be used in the component combinations approved by the manufacturer (see instructions for use and catalogues). The manufacturer is not liable for damage caused by component combinations and applications that have not been approved by the manufacturer. Opening and repairing of the products may only be carried out by Macu or with the written permission of Macu.

47. Conformity

The Explorer socket complies with the general safety and performance requirements set out in Regulation (EU)2017/745 for medical devices.

Legal manufacturer

macu4 AG, Rämistrasse 18, 8001 Zurich, Switzerland, CHE-301.889.074

Special production

Macu declares that they comply with the relevant European standards for the design, manufacture and supply of custom prosthetic products. Continuous compliance with the standard is monitored through a programme of internal and external audits.

The socket of the Explorer System meets the requirements of the European Regulation (EU)2017/745 Annex I for medical devices. Based on the classification criteria according to Annex VIII and Annex XIII of this regulation, the product has been classified as Class I and custom-made device. The declaration of conformity has therefore been drawn up by the manufacturer under its sole responsibility in accordance with Annex IV to the Directive.

Trademark


All designations mentioned in this document are subject without restriction to the provisions of the applicable trademark law and the rights of the respective owners.

All trademarks, trade names or company names mentioned herein may be registered trademarks and are subject to the rights of their respective owners.

The absence of explicit identification of the trademarks used in this document does not imply that a designation is free of third-party rights.

48. Symbols

The table below shows all symbols (according to ISO 15223-1:2021) used on the product and on the label.

Custom made	
	This symbol indicates that a manufacturer declares that the product meets all legal requirements for custom-made production.
Consult operating instructions or electronic instructions for use	
	This symbol indicates that the user must consult the instructions for use.
Manufacturer	
	This symbol indicates the manufacturer of the medical device.

Warnings	
	This symbol indicates that caution is required when operating the product or that the current situation requires the operator's attention or action to avoid undesirable consequences.
Temperature limit	
	This symbol indicates the temperature limits to which the medical device can be safely exposed.
Date of manufacture	
	This mark indicates when the medical device was manufactured.
Batch designation	
	This mark indicates the manufacturer's batch code so that the batch or lot can be identified
Catalog number	
	This mark indicates the manufacturer's catalog number so that the item can be identified.
Medical device	
	Indicates that the item is a medical device.
Keep away from sunlight	
	This symbol indicates a medical device that must be protected from light sources.
Magnetic fields	
	This symbol indicates that magnets are contained in the product. The magnetic field of magnets may affect or disable the function of medical devices (e.g. pacemakers or defibrillators).

49. Technical data

49.1. Socket

General		
Drawing	Relatively small socket:	Relatively large socket:
		
Dimensions	n/a	
Net weight (without lacing system, screw nut, socket ring and padding materials)	Between approx. 0.028 kg / 0.062 lbs (relatively small custom-made stock) and 0.1 kg / 0.22 lbs (relatively large custom-made stock)	
Lifetime of the 3D printed socket base component	2 years at an average of 120 minutes per day of daily use on the 3D printed socket component supplied by Macu.	

Maximum load forces in the coupling interface	
Axial tension and compression force	ca. 250 N
Bending moment around socket axis	ca. 5 Nm
Bending moment normal to socket axis	ca. 10 Nm

Environmental conditions	
Storage and transport in the original packaging	0 °C/32 °F to +35 °C/95 °F max. 70% relative humidity
Storage and transport without original packaging	0 °C/32 °F to +35 °C/95 °F max. 70% relative humidity
Operation	0 °C/32 °F to +35 °C/95 °F max. 100% relative humidity
Storage	0 °C/32 °F to +35 °C/95 °F max. 70% relative humidity and protect from direct sunlight.

Explorer Socket		
01011	Explorer Socket	Connector M12x1.5'
01012	Explorer Socket	Connector 1/2"-20

50. Returns

If you believe there is a problem with your Explorer Socket, please first read this manual carefully before contacting support@macu4.com or your prosthesis manufacturer (unless you purchased the product directly from Macu). Please use the Product Experience Report (PER): https://25462115.fs1.hubspotusercontent-eu1.net/hubfs/25462115/mac4_AG_Explorer/Product_Experience_Report_EN.pdf

Please include the item and lot number (one single number) of your stock when requesting warranty repairs or returns. You will find these printed on the inside of the stock. If you are returning the socket to us, please pack it appropriately; the original packaging in which you received your module is excellent for this purpose.

51. Disposal of products

If you are no longer going to use a product, you can dispose of it as follows:

- Return to macu4 AG (for clients from CH)
- See FAQs on the website Website www.macu4.ch

52. Document history

Version	Date	Description
1.0	28.12.2021	Creation
2.0	07.02.2022	Update
3.0	01.05.2022	Update
4.0	-	-
5.0	30.12.2022	Update Chepter 9 and Chapter 13
6.0	15.04.2023	Update Chepter 9
7.0	22.08.2023	Update Chapter 2
8.0	01.02.2024	Update Chapter 7, 9, 12, 13