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Instruction for Use | EN

1. Vorwort

Gültig ab: 25.08.2025

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1. Foreword

Valid from: 25.08.2025

- Read this document carefully before using the product.
- Follow carefully the warnings and instructions for assembly to avoid product damage.
- Medical personnel should assemble the product.
- Contact the manufacturer if you have any questions about the product (e.g., commissioning, use, maintenance, unexpected operation, or incidents).
- Keep this document in a safe place.

These instructions for assembly provide you with important information on the use, adjustment, and handling of the provided component(s).

2. Product

Product service description

The Lynk Interface is a semi-finished part which can be integrated into a custom-made orthosis or prosthetic socket. It must be only used together with Lynk Modules.

Our role is limited to supplying the Interface as a non-finished component. The medical personnel is solely responsible for the integration of this component with other materials or parts, the assembly of the complete prosthetic device, and its clinical evaluation and fitting for the end user. In this capacity, the medical personnel acts as the Legal Manufacturer of the final medical device under applicable medical device regulations and assumes full responsibility for regulatory compliance, safety, and performance of the finished device.

CAUTION

3. Warnings

FAILURE TO HEED THESE WARNINGS OR TO FOLLOW THE INSTRUCTIONS FOR ASSEMBLY COULD RESULT IN MALFUNCTION OF THE PRODUCT OR PHYSICAL INJURY. IN THE EVENT OF INJURY, SEEK MEDICAL ATTENTION PROMPTLY.

Notes on unpacking the supplied component

Only open the packaging at the indicated locations with care so that the component inside is not damaged. Inspect the delivered component and be sure that it is undamaged.

Notes on the lamination process

Make sure that the Interface does not deform, for example, when you work with high temperatures or high pressures applied to the Interface.

Notes on using the connector

Make sure that you do not weaken the Interface, for example, when you decide to use cutting or abrasive processes. Make sure that no material remains inside the Interface before using it together with a Lynk Module.

4. Approach

Material needed


Interface	1x assembly material set (content: Interface for screwing, 8x Screw (M4x20), 8x Square nut (M4))
Additional materials needed on behalf of medical personnel	Allen wrench

Example to illustrate how to use the component

Phase 1: Planning & Preparation


Determine the Interface Position

- Use the user's arm or a plaster cast to simulate the final shape of the orthosis.
- Attach a module to the Interface to visualise its position on the orthosis.
- Ensure sufficient clearance for movement and additional orthosis features (e.g., wrist immobilisation).
- Mark the desired location on your design before fabrication.

Design the Orthosis with Mounting Considerations


- Ensure a flat surface to mount the Interface.
- Decide how the M4 square nuts will be integrated:
 - Option 1: Embed the nuts in the orthosis material during fabrication.
 - Option 2: Drill holes and attach the nuts after fabrication (e.g., glueing or securing in a housing).
- Plan how to protect the user's skin from contact with the nuts and screws (e.g., padding, embedded housing).

Phase 2: Orthosis Fabrication


Create the Cuff Structure

- Fabricate the cuff according to the planned design, ensuring structural integrity and user comfort.
- If embedding nuts during fabrication, securely position them within the material at the designated locations. Note that in this case, the holes for the screws shall also be considered to prevent damaging the nuts afterwards.
- If the nuts will be added later, ensure there is enough space for proper fixation without compromising the cuff's strength.

Mark & Drill the Screw Holes


- Use the Interface as a guide or follow the provided dimensions to mark the hole positions.
- Drill holes carefully, ensuring they align precisely with the embedded or future nut placements.
- If the nuts were not embedded earlier, insert and secure them now using the preferred method (e.g., glueing, press-fitting).
- Check that the holes and nuts allow for smooth screw insertion without misalignment.

Phase 3: Assembly & Testing


Attach the Interface to the Orthosis

- Position the Interface on the pre-drilled holes.
- Secure with 8x M4x20 screws or adjust the screw length as needed based on your orthosis's design.
- Confirm that the nuts do not protrude uncomfortably or affect the orthosis's structure.

Final Check & Functionality Test

- Insert a module into the Interface and check for interference with the orthosis.
- Verify that the Interface is securely attached and can handle intended activities.
- Ensure no hardware (nuts, screws) is exposed to the user's skin.
- Adjust or add cushioning if necessary.

Note: The illustrations provided are for indicative purposes only and may not be to scale.

5. Technical data

General	
Drawing	<p>Proximal end Distal end</p>
Dimensions	<p>Dimensions: 110, 105, 71.67, 38.33, 5, 51.4, Ø4 X8, R4 X4, 11.59, 2.5, Ø8 X8, A-A</p>
Net weight	Between approx. 30 g / 1.05 oz
Material of printed component	Module body: PA2200, coloured with DM standard colour (DyeMansion)
Lifetime of the 3D printed component	2 years at an average of 120 minutes per day of daily use on the 3D printed component
Material of the Nut and screw	Stainless steel A2
Dimensions of the threading	Metric: M4x20'

6. Legal notice

Limitation or exclusion of liability

Failure to follow these instructions for use, care, maintenance, handling, storage, etc., of this product, or the use of this product in any component combination other than those approved by macu4 AG (hereinafter referred to as “Macu”) constitutes misuse of the product and invalidates any and all warranties, express or implied, in their entirety as to any consequences of such misuse. In the event of such misuse, Macu disclaims liability for any adverse consequences to the maximum extent allowed by law. Any opening, disassembly, or repair of the products may only be carried out by Macu or with Macu’s written permission.

Manufacturer of the semi-finished component

Macu4 AG, Rämistrasse 18, 8001 Zurich, Switzerland

Trademark

All designations mentioned in this document are subject without restriction to the provisions of the applicable trademark law and the rights of the respective owners. All trademarks, trade names or company names mentioned herein may be registered trademarks and are subject to the rights of their respective owners. The absence of explicit identification of the trademarks used in this document does not imply that a designation is free of third-party rights.

7. Environmental conditions

Conditions	
During transport	<ul style="list-style-type: none">• Protect from direct sunlight (UV radiation)• Temperature shall be 0 °C to +35 °C / 32 °F to 95 °F

8. Allergic reactions

THIS PRODUCT IS FORMULATED FOR GENERAL SKIN COMPATIBILITY. HOWEVER, ALLERGIC REACTIONS MAY OCCUR IN SOME INDIVIDUALS, PARTICULARLY THOSE WITH A HISTORY OF SKIN ALLERGIES OR SENSITIVITY. THE PRODUCT COMPONENTS MUST NOT BE USED IF THE INDIVIDUAL HAS A KNOWN ALLERGY OR HAS PREVIOUSLY EXPERIENCED A REACTION TO ANY OF THE MATERIALS USED IN THIS PRODUCT. DISCONTINUE USE IMMEDIATELY AND SEEK MEDICAL ADVICE IF REDNESS, ITCHING, RASH, SWELLING, OR IRRITATION OCCURS. INDIVIDUALS WITH SENSITIVE SKIN SHOULD REVIEW THE FULL LIST OF MATERIALS AND, IF NECESSARY, PERFORM A PATCH TEST BEFORE FIRST USE. FOR DETAILED INFORMATION ON THE MATERIALS USED IN SPECIFIC PRODUCTS, PLEASE CONTACT MACU. THE TECHNICAL DATA SECTION PROVIDES GENERAL INFORMATION ON MATERIAL USAGE ONLY.

9. Returns

If you suspect or believe there is a problem with the semi-finished component, please first read this manual carefully before contacting support@macu4.com.. Please use the Product Experience Report (PER):

https://25462115.fs1.hubspotusercontent-eu1.net/hubfs/25462115/mac4_AG_Lynk/Product_Experience_Report_EN.pdf