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1. Instructions for Use (EN) | Lynk System

Valid from: 01.02.2024

- Read this document carefully before using the product.
- Observe the safety instructions to avoid injuries and product damage.
- Instruct the user on the proper and safe use of the product.
- Contact the manufacturer if you have any questions about the product (e.g., commissioning, use, maintenance, unexpected operation or incidents). You will find the contact details on the back.
- Keep this document in a safe place.

These instructions for use provide you with important information on the use, adjustment and handling of the Lynk Ortho-prosthetic System. **Please notice that in addition to this Instructions for use, there are several explanation videos available. See <https://macu4.com/en/application-lynk>.**

1. Use

The Lynk system can generally be used when a person wants to support or protect the hand function and/or wrist from loading forces during an activity that is based on a gripping function with, for example, pulling or pushing forces.

Intended Use


The Lynk System is intended for persons with partial or complete loss of function of the hand due to temporary or permanent weakened hand or hand injury, paralysis, partial hand absence or wrist disarticulation absence. By using the Lynk Cuff with the appropriate accessory, a lack of gripping dexterity, capability or force can be replaced or supported.

Indications

- The system is suitable for persons with impaired or missing gripping function.
- The Ortho-Prosthetic system can be used for people with unilateral congenital hand deficiencies or amputation immediately above the wrist (according to ISO 8548-2:2020) and, depending on the arm stump shape.
- The system is suitable for persons two years and older. An individual evaluation by the specialist is recommended to assess whether the user is physically and mentally in the condition to use the system properly.

Contraindications

- The Ortho-Prosthetic System is not indicated for individuals two years and younger. Exception: The specialist (e.g. doctor, therapist, orthopaedic technician) can make a proper assessment of whether the user is in a condition to use the system anyway.
- The Ortho-Prosthetic System is not indicated for individuals with bilateral forearm deficits. Exception: The specialist (e.g. doctor, therapist, orthopaedic technician) can make a proper assessment of whether the user is in a condition to use the system partially anyway.
- The Lynk System is not indicated for individuals with a severe limitation of the upper limb. Exception: The specialist (e.g. doctor, therapist, orthopaedic technician) can make a proper assessment of whether the user is in a condition to use the system partially anyway.
- The Lynk Cuff is not indicated for individuals with fresh wounds, allergies or sensitive skin on the arm.
- The Lynk Cuff is not indicated for individuals with a circumference of the forearm larger than 330mm, a circumference of the wrist smaller than 130mm, and a forearm length shorter than 140mm.

- The system may only be used by users who are physically and mentally fit to operate the system properly. If you are unsure whether the physical and/or mental fitness is sufficient, we recommend an individual assessment by a specialist (e.g. doctor, therapist, prosthetist).

Lynk System

Combination possibilities

The system has been tested for combination possibilities of the Lynk Cuff with the different Lynk modules. When using the Lynk Cuff and to achieve an improved hold on the forearm by an additional upper arm connection, the professional has to check on his own responsibility if the combination possibility is given. When using the Lynk Interface as a semi-finished part for a self-manufactured orthosis or cuff, the specialist must check on his/her own responsibility whether it can be combined. For a cuff other than ours, the specialist must check on his/her own responsibility whether it can be combined with the Lynk modules.

Operating conditions - General

The Lynk System was developed especially for everyday, sports and leisure activities with low to moderate stress. The system must not be used for extraordinary activities. These extraordinary activities include, for example, sports with excessive coupling and/or impact loads (push-ups, downhill and mountain biking, etc.) or extreme sports (free climbing, paragliding, etc.). In addition, the system should not be used for driving motor vehicles, driving heavy equipment (e.g., construction machinery), operating industrial machinery or operating motor-driven work equipment.

Operating conditions – Push Module

To achieve the highest safety level, the Push Module was tested mechanically based on the *EN ISO 11334-1:2007 Assistive products for walking manipulated by one arm - Requirements and test methods - Part 1: Elbow crutches*. According to the requirements of the above mentioned ISO, the use of the Push Module to support the walking with crutches is limited to a maximum body weight of 65kg. A higher body weight in combination of the above mentioned operating condition will potentially lead to material deformation, inducing instability at the interface between the module of the crutch. Specialists should contact support@macu4.com for further information.

**2. Warnings**General

Carefully review your components before using them to make sure that no area of the component is not damaged. Please always contact support@macu4.com and don't use the component(s) when you have the feeling that, for instance

- a component is damaged.
- a component may not be suitable for an indicated activity.
- a module cannot be inserted easily into the Cuff insertion slot.
- the Velcro or Lacing closure is not working properly or releases unintentionally.

Make sure to watch the instructional video for your intended activity. Incorrect use may cause inconvenience or an accident. If you don't find your activity among the instructional videos provided, please contact support@macu4.com or the professional who provided you with the Lynk System for instructions before starting the activity.

Notes on Modules - General

Make sure that you do not overload the Modules. Make sure that you always clean the Modules after use. Remove visible dirt particles so that Modules with flexible modes are not restricted in their function.

Don't use any Module as a means to manage objects that are sharp or that could cause damage in case of loss of control (i.e. raser, drill, etc.).

Don't use a Module if the diameter of the object is too large (see the specifications for each module). If the diameter of an object is too small, you can increase the diameter slightly by using a soft material around the area where you want to attach the module to the object. For example, when using the Hold Module to hold a fork.

Before using a Module, check if it is compatible with the handle (i.e. crutch handle, bike handle, etc.) or object (i.e. fork, glass, etc.), that the Module is tightly secured on the handle or the grip or the object, that you have chosen the correct orientation of the Module on the object/handle.

Notes on CLIC Module

Make sure that the Clic Head is correctly placed in the Clic Adapter and that you tested prior to usage how to connect and disconnect either the module from a handle or the Clic Adapter from the Clic Head. Before using any device with the Clic Module, make sure that you can properly connect and disconnect the adapter to and from the Clic Head. The Clic Module is not suitable when you need to pull on the Module while performing your activity (i.e. jumping with a scooter) because the Clic Adapter disconnects from the Clic Head when you apply pull forces.

Notes on HOLD Module

Make sure that the lever of the ball joint is tightened and that the Velcro strap is correctly tightened to hold the object well enough and not come loose unintentionally. Do not overtighten the eccentric cam, as this may damage the material.

Notes on HOOK Module

Make sure that the handle is correctly placed in the Module and that you respect the maximum allowed load cases.

Notes on PUSH Module

Make sure that the Module is correctly connected to the handle of the object (e.g. crutch) and that you respect the maximum allowed load cases. See additionally the Operating Conditions "Push Module".

Before starting any activity, make sure that the module head (black curved handle element) is correctly attached to its module base (it should not be rotatable). If this is not the case, the screw may have loosened slightly. In this case, tighten the screw.

Notes on selection of product size

Don't use a Cuff size that is too big for your arm. Any excessive tightening can lead to a dysfunction of the Cuff Hook(s).

Don't use a Cuff size that is too small for your arm. It can lead to painful pressure marks and/or impede vascular flow.

Notes on using the Cuff and Cover

Don't use the Cuff with/without Cover when the protective material is missing.

Make sure that the Velcro straps respectively the Lacing Knob (when using the Cover) are facing upside (posterior) on your forearm.

Make sure that the Cuff (with or without Cover) always sits tight on your forearm. If it gets loose, tighten the Velcro respectively Lacing system.

Notes on ambient conditions

Make sure that you store or use the products in accordance with the information on temperature, UV radiation and humidity given in Chapters 9 and 10. The material properties of the components may change significantly when used and stored at higher or lower temperatures, as stated in this document, and lead to inconvenience or an accident.

Notes on storage conditions

Make sure that the components are well protected during transportation and permanent storage. Protect especially the Button area of each Module from load situations (e.g. don't place heavy objects on it or apply bending forces). It could deform and/or damage the Module interface. In case of deformation or damage, the Module can no longer be used with the cuff.

3. Use of the Components

Before using the Lynk components for the first time, familiarise yourself with the instructions for use. The following chapters explain how to use the Lynk Modules together with the Lynk Cuff.

3.1. How to use the Lynk Cuff

Design and structure of the Lynk Cuff

Cuff from the top side equipped with Velcro System

Cuff from the backside (here: with two Hooks in their anchoring position)

Cuff from the top side equipped with Cuff Cover and Lacing System

3.1.1. How to use the Lynk Cuff with Velcro System

When using the Cuff for the first time

Make sure that the Cuff sits tight and cannot be moved easily. If the Cuff sits too loose on your arm, you should adapt the Velcro straps. The fit of the Cuff on the arm should be snug but not uncomfortable.

Option 1 | Open the Velcro strap at the level of the Velcro part that shows the macu4 logo. Don't remove the part that shows the logo. Pull slightly on the Velcro strap to tighten the sit on your arm.

Option 2 | Open the Velcro strap at the level of the Velcro part that shows the macu4 logo. Remove the part that shows the logo. Shorten the Velcro strap with scissors. Reattach the macu4 logo Velcro part. Pull slightly on the Velcro strap to tighten the sit on your arm.

Put on the Lynk Cuff

Insert your arm | You can either place the Cuff on a table and insert your arm into it. Or you can use your other hand to hold the Cuff while inserting your arm.

Secure the Cuff on your arm | Place the Hook in front of the Hook Anchoring Holes - the magnets will allow it to click into place. Make sure it is well positioned and secured by gently moving the Cuff.

3.1.2. How to use the Lynk Cuff with Lacing System

When using the Cuff Cover

Whenever you want to use the Cuff together with the Cuff Cover, you must remove the Velcro System. Store the Velcro System in a safe place and reinstall it to the Cuff whenever you need it.

Release the Velcro closure system | The Velcro fastener is attached to the Cuff on one side via slits. In this area, you can release the Velcro closure.

Remove the Velcro closure system | Carefully pull the Velcro straps out of the slits. Once released, store the Velcro closure system in a secure place to use it at a later time again, if needed.

Put on the Lynk Cuff

Note: Whenever you want to use the Cuff Cover, the Velcro closure System has to be removed, as explained above.

Insert your arm | You can either place the Cuff on a table and insert your arm into it. Or you can use your other hand to hold the Cuff while inserting your arm.

Prepare Cuff Cover | Take the Cuff Cover and make sure that the Lacing Knob is in the released position by softly pulling it up.

Anchor the Hooks | Place the Hook in front of the Hook Anchoring Holes. The Hook is equipped with magnets, and this will allow it to click into place.

Tighten the Cuff | When both Hooks are in place, press the Lacing Knob and turn it clockwise. Turn until the Cuff fits snugly on the arm. The fit of the Cuff on the arm should be snug but not uncomfortable.

3.1.3. Final control

Check the fit of the Cuff before you start any activity. Make sure it is well positioned and secured by trying to move the Cuff gently. The Cuff should remain firmly in position. If the Cuff is not tight and can be moved, recheck the Lacing System or Velcro System.

3.1.4. How to remove the Lynk Cuff

Option 1 | Applicable for both closure systems. Pull on the Hook. It shall disengage from the Cuff. If the Hook is free, you can remove the Cuff from your arm.

Option 2 | Applicable for lacing system only. Alternatively, you can first loosen the Lacing System by pulling slightly on the Knob before pulling on the Hook.

3.2. How to use a Lynk Module with the Lynk Cuff

Insert Module | Insert the module into the interface section through the slit in the cuff as shown. The button should face you so that it remains visible in the interface section. You should hear a small sound.

Change position | Press the Button of the Module to change its position in the Cuff and to adapt it to your needs and comfort. Press the Button as well when you want to remove the Module.

Final check | Make sure the Module is well positioned and secured by trying to move it gently. The Module should remain firmly in position.

3.3. CLIC Module

Design and structure of the Lynk CLIC Module

The Clic Module is provided as a set of 1 Clic Adapter, assembled with 3 different Clic Lamellas and 3 different Clic Heads.

Disassembled | Module with Clic Head uncoupled from Clic Adapter (Clic Lamella is inserted into Adapter)

Assembled | Module with Clic Head coupled to Clic Adapter

Lynk Clic Module Heads

The Clic Module is offered with different Clic Heads to support the user in a range of activities that may require different rotation levels of the Clic Head.

- **Clic Head 360** | Module Head providing a complete rotation, offering complete freedom.
- **Clic Head 70** | Module Head allowing a partial rotation to prevent constraining the arm movement completely.
- **Clic Head 0** | Module Head offering no rotation but various locked positions. You can choose between different locked angles.

Lynk Clic Lamellas

The Clic Module is offered with 3 different Clic Lamellas to support the user in a range of activities that may require different release forces when attaching and detaching a Clic Head. The Clic Module is assembled with Clic Lamella 1 by default.

- **Clic Lamella 1** | Lamella for lower forces needed to connect and disconnect Clic Head and Clic Adapter.
- **Clic Lamella 2** | Lamella for moderate forces needed to connect and disconnect Clic Head and Clic Adapter.
- **Clic Lamella 3** | Lamella for higher forces needed to connect and disconnect Clic Head and Clic Adapter.


You can change the Lamella type as follows:

Remove the Pin | Press on the pin as illustrated with a small tool (no scarf tip) such as a pen onto the pin and in the direction of the Modu until the Pin is released.

Change the Lamella | Slide out the Clic Lamella on the side through the slit. Insert the new Lamella. Make sure that the pin hole remains free.

Insert the Pin | Insert the pin to secure the new Clic Lamella. If you cannot insert it with your hand, gently press a flat object (i.e. book) against the module.

3.3.1. Connect the Clic Module to a handle | Option 1

Connect the Clic Module to a handle by pushing the clamp area of the Clic Head onto the bar. Release it by pulling. Check the fit of the Clic Module on the handlebar. If the Clic Head is too loose on the handlebar, you can add some rubber material to the handle to improve the grip. The supported diameter range of a handle is provided in Chapter 10. Always test to connect and disconnect the Clic Module to and from the handle before starting any activity.

- **Connect to handle** | Connect the Clic Module to a handle by pushing the clamp area of the Clic Head onto the bar. Find the most suitable location on the handle.
- **Disconnect from handle** | Disconnect the Clic Module from the handle by pulling it from the handle.

3.3.2. Connect the Clic Module to a handle | Option 2

Instead of attaching and detaching the entire module to and from the handle, you can attach the Clic Head permanently to the handle by using cable ties, as explained below. The supported diameter range of a handle is provided in Chapter 10. To use the Clic Module in this option 2, you connect to the Clic Head via the Clic Adapter. Insert the Clic Adapter onto the Clic Head. And releasing it by pulling. Always test to connect and disconnect the Adapter to and from the Clic Head before starting any activity.

- **Attach Clic Head** | Connect the Clic Head to a handle by pushing the clamp area of the Clic Head onto the bar. Find the most suitable location on the handle.
- **Secure Clic Head** | If the Clip Head sits loose, add some rubber material to the handle to improve the grip. Insert a Velcro strap or cable ties through the slit of the Clip Head and tighten them.
- **Connect Clic Adapter & Head** | Connect to the Clic Head that is fixated to the handle by pushing the Adapter slit onto the Clic Head. And releasing it by pulling on the Clic Head.

3.3.3. Use case examples with the Module

Please consult the website of macu4, where you will find specific explanation videos under the section 'Application'.

- **Using a broom** | Attach the Clic Module to a broom. Test yourself which Clic Head works best for you.
- **Using a bike** | Attach the Clic Module to the handlebar of a bike. Test yourself which Clic Head works best for you.
- **Using a scooter** | Attach the Clic Module to a broom. Test yourself which Clic Head works best for you.

3.4. HOOK Module

Design and structure of the Lynk Hook Module

The Hook Module is provided as one single part.

View from front | The Hook Module has a slim design while made from a stable structure. The design of the hook part includes the opening towards the cuff.

View from side | The shape of the curved hook section allows objects with weight to be pulled. The overlapping plateau element allows additional weight to be pushed.

Apply pushing force | The overlapping plateau element at the other end of the curved hook section allows additional weight to be pushed.

Use case examples with the Module

Please consult the website of macu4, where you will find specific explanation videos under the section 'Application'. The maximum handle diameter and object weight are provided in Chapter 10.

Holding a bag | Insert the handle of the bag into the Hook.

Pulling a suitcase | Insert the handle of the suitcase in the Hook and start pulling gently.

Pushing and pulling weight in the gym | Clip your module to the Handle of the fitness equipment and simply push it. Release the Module from the bar by turning it while pulling on it softly.

3.5. Push Module

Design and structure of the Push Module

The Push Module is provided as one single part with a removable Velcro element.

Flexible module head | The module head is a curved grip element and is made from a flexible material to better conform to a handle.

Velcro closure | A Velcro strap is connected to the module via a slit in the head and the Velcro hook. The hook is held in position thanks to the magnets and a click-in mechanism.

Adaptable orientation | You can slightly turn the module head to the left and right to find the desired orientation for your handle.

Use case examples with the Module

Always make sure that objects or handles are properly positioned and secured by lifting the object and shaking it gently. Please consult the website of macu4, where you will find specific explanation videos under the section 'Application'.

- **Using a crutch** | Make sure the Velcro strap is free and place the grip part onto the handle. Turn the Velcro strap around the handle and secure the Velcro hook in its position. If needed, shorten the length of the Velcro strap.
- **Attaching cylindrical objects** | Make sure the Velcro strap is free. Press the object onto the grip part. Turn the Velcro strap around the object and secure the Velcro hook in its position. If needed, shorten the length of the Velcro strap or replace the Velcro strap with a longer one.
- **Attaching to other objects** | Use the Push module to fixate an arch and perform arch shooting. Proceed as described for the other used cases.

3.6. Hold Module

Design and structure of the Hold Module

The Hold Module is provided as one single part with a removable module head element.

Flexible module head | The module head is a clamp-like element and is made from a flexible material to better conform to objects.

Velcro closure | A Velcro strap is attached to the module head through a slot in the head. The strap can be wrapped around the module head and fixed in place.

Adaptable orientation | You can turn the module head in various directions to find the desired orientation for your activity.

Use case examples with the Module

Always make sure that objects or handles are properly positioned and secured by lifting the object and shaking it gently. Please consult the website of macu4, where you will find specific explanation videos under the section 'Application'.

- **Kitchen accessories** | Position the kitchen aid object into the Hold Module Head. Secure its position by tightening the Velcro strap.
- **Flat objects** | Insert the flat object into the lateral slit that you find on the Hold Module Head.

4. Maintenance

Dirt particles in the spaces between the products should be removed with a soft-bristle brush. This pre-cleaning of the products should be done without the use of water.

Product surfaces should be cleaned with a damp cloth (non-scratching). Cold or lukewarm water should be used for this purpose.

5. Adverse event reporting

Any serious incident related to the device should be reported to the manufacturer and the competent authority of the member state where the user and/or patient is established. Please use the form Please use the Product Experience Report (PER)

https://25462115.fs1.hubspotusercontent-eu1.net/hubfs/25462115/mac4_GmbH_Lynk/Product_Experience_Report_EN.pdf and send it to by email to support@macu4.com.

6. Legal notice

Limitation or exclusion of liability

macu4 GmbH, hereinafter referred to as Macu, is only liable if the specified handling and processing instructions, as well as the care instructions and maintenance intervals for the product, are observed. Macu expressly points out that this product may only be used in the component combinations approved by the manufacturer (see instructions for use and catalogues). The manufacturer is not liable for damage caused by component combinations and applications that have not been approved by the manufacturer. Opening and repairing of the products may only be carried out by Macu or with the written permission of Macu.

7. Conformity

The product components of the system Lynk with the Basic UDI-DI marking

Lynk Cuff	426242965P020201R4
Lynk Cover	426242965P020202R6
Lynk Clic	426242965P020207RG
Lynk Hold	426242965P020206RE
Lynk Hook	426242965P020204RA
Lynk Push	426242965P020205RC

meet the general safety and performance requirements according to Regulation (EU)2017/745 for medical devices and the Medical Devices Ordinance (MedDO) of 1 July 2020 (812.213).

Legal manufacturer

macu4 GmbH, Bücklestrasse 3, 78467 Constance, Germany

CE Marking Details

Macu declares that it complies with the relevant European standards for the design, manufacture and supply of ortho-prosthetic products. The Modules of the Lynk System meet the requirements of the European Regulation (EU)2017/745 Annex I for medical devices. Based on the classification criteria according to Annex VIII of this regulation, the product has been classified as Class I “Accessories”. The declaration of conformity has, therefore, been drawn up by the manufacturer under his sole responsibility in accordance with Annex IV of the Regulation. All individual Modules of the Lynk System are marked accordingly.

Trademark

All designations mentioned in this document are subject without restriction to the provisions of the applicable trademark law and the rights of the respective owners. All trademarks, trade names, or company names mentioned herein may be registered trademarks and are subject to the rights of their respective owners. The absence of explicit identification of the trademarks used in this document does not imply that a designation is free of third-party rights.

8. Symbols

The table below shows all symbols (according to ISO 15223-1:2021) used on the product and on the label.

CE Mark

This symbol indicates that a manufacturer declares that the product meets all the legal requirements for CE marking and can be sold throughout the EEA.

Instructions for use.

This symbol indicates that the user must consult the instructions for use.

Manufacturer

This symbol indicates the manufacturer of the medical device.

Warnings

This symbol indicates that caution is required when operating the product or that the current situation requires the operator's attention or action to avoid undesirable consequences.

Temperature limit

This symbol indicates the temperature limits to which the medical device can be safely exposed.

Date of manufacture

This symbol indicates when the medical device was manufactured.

Batch designation

This symbol indicates the manufacturer's batch code so that the batch or lot can be identified.

Catalog number

This symbol indicates the manufacturer's catalogue number so that the item can be identified.

Medical device

This symbol indicates that the item is a medical device.

Clear device identification

This symbol identifies a carrier that contains information for unique device identification.

Keep away from sunlight

This symbol indicates a medical device that must be protected from light sources.

Magnetic fields


This symbol indicates that magnets are contained in the product. The magnetic field of magnets may affect or disable the function of medical devices (e.g. pacemakers or defibrillators).

9. Technical data


9.1. CUFF

General	
Drawing	
Dimensions (size: short)	Small: 83mm x 61mm x 127mm Medium: 93mm x 68.5mm x 130.5mm Large: 105mm x 76.5mm x 132mm
Dimensions (size: standard)	Small: 87.5mm x 66mm x 158mm Medium: 97mm x 72mm x 158mm Large: 112mm x 81mm x 158mm
Weight (size: short, assembled)	Small: 67 g Medium: 70 g Large: 75 g
Weight (size: standard, assembled)	Small: 75 g Medium: 80 g Large: 87 g
Material	Cuff shell: PA2200, coloured with DM standard colour (DyeMansion) Protective material: Alcantara, Velour Velcro: Nylon Magnets: Neodymium
Lifetime of the component	2 years at an average of 120 minutes per day with daily use

9.2. COVER

General	
Drawing	
Dimensions	Small: 83.5mm x 39.5mm x 131mm Medium: 93.5mm x 44.5mm x 132.5mm Large: 105mm x 50mm x 134mm
Weight (assembled)	Small: 37 g Medium: 40 g Large: 43 g
Material	Cuff shell: PA2200, coloured with DM standard colour (DyeMansion) Protective material: Alcantara, Velour Lacing System: Fitgo Magnets: Neodymium
Lifetime of the component	2 years at an average of 120 minutes per day with daily use


9.3. CLIC Module

General	
Drawing (Clic Module)	
Dimensions (Clic Module)	167mm x 45.5mm x 44mm
Weight (Clic Module)	57 g
Drawing (Heads only)	
	Head 360°
Drawing (Heads only)	
	Head 70°
Drawing (Heads only)	
	Head 0°
Dimensions (Head)	44mm x 44mm x 44mm
Weight (Head)	11.5 g
Drawing (Lamella)	
	Lamella soft
Drawing (Lamella)	
	Lamella medium
Drawing (Lamella)	
	Lamella strong
Dimensions (Lamella only)	39mm x 35.3mm x 3.4mm
Drawing (Adapter only)	
Dimensions (Adapter only)	44mm x 167mm x 10mm
Weight (Adapter)	23 g
Allowed handle diameters	From 26 mm to 32 mm

Material	Module body: PA2200, coloured with DM standard colour (DyeMansion)
Lifetime of the Module	2 years at an average of 120 minutes per day with daily use

Maximum load forces in the Module	
Axial pushing force	Approximately 200 N

9.4. HOLD Module

General	
Drawing	
Dimensions	39.5mm x 198mm x 59mm * *Dimensions may vary depending on the orientation of the Module
Allowed diameters for cylindrical objects	From 10mm to 85 mm * *Dimensions may vary depending on the shape and material of the object.
Allowed thickness of flat objects	From 0.3mm to 12mm * *Dimensions may vary depending on the shape and material of the object.
Weight	86 g
Material	Module body: PA2200, coloured with DM standard colour (DyeMansion) Module head: Ultrasint TPU01
Lifetime of the Module	2 years at an average of 120 minutes per day / daily use


Maximum conditions the Module can withstand	
Maximum object weight	3 kg
Maximum pushing force	30 N
Maximum pulling force	30 N

9.5. HOOK Module

General	
Drawing	
Dimensions	49.2mm x 168.8mm x 29.1mm
Allowed diameters for handles	Up to 32 mm
Weight	30 g
Material	PA2200, coloured with DM standard colour (DyeMansion)
Lifetime of the Module	2 years at an average of 120 minutes per day / daily use

Maximum conditions the Module can withstand	
Maximum object weight	50 kg
Maximum pushing force	200 N
Maximum pulling force	400 N

9.6. PUSH Module

General	
Drawing	
Dimensions	49.2mm x 168.8mm x 29.1mm
Allowed diameters for handles	From 30mm to 40mm * *Dimensions may vary depending on the shape and material of the object.
Weight	58 g
Material	Module body: PA2200, coloured with DM standard colour (DyeMansion) Module head: Ultrasint TPU01
Lifetime of the Module	2 years at an average of 120 minutes per day / daily use

Maximum conditions the Module can withstand	
Maximum body mass when using crutches	65 kg
Maximum pushing force	500 N
Maximum pulling force	100 N

9.7. Environmental conditions

Environmental conditions	
Storage and transport in the original packaging	0 °C/32 °F to +35 °C/95 °F max. 70 % relative humidity
Storage and transport without original packaging	0 °C/32 °F to +35 °C/95 °F max. 70 % relative humidity
Operation	0 °C/32 °F to +35 °C/95 °F max. 100 % relative humidity
Storage	0 °C/32 °F to +35 °C/95 °F max. 70 % relative humidity and protect from direct sunlight.

10. Configuration of the System

Please see the measuring instructions on the macu4 website under the heading Lynk Modules.

11. Returns

If you believe there is a problem with your product, please first read this manual carefully before contacting support@macu4.com or your prosthesis manufacturer (unless you purchased the product directly from Macu). Please use the Product Experience Report (PER): https://25462115.fs1.hubspotusercontent-eu1.net/hubfs/25462115/mac4_GmbH_Lynk/Product_Experience_Report_EN.pdf

Please provide the UDI-DI (starting with “(01) ____”) and UDI-PI (starting with “(10) ____”) of your Module when requesting warranty repairs or returns. These can be found imprinted on the Module.

If you return your product to us or your prosthetist, please pack it appropriately; the original packaging in which you received your Module is perfectly suitable for this purpose.

12. Disposal of products

If you are no longer going to use a product, you can dispose of it as follows:

- Return to macu4 GmbH
- See FAQs on the website Website www.macu4.com

13. Document history

Version	Date	Description
---------	------	-------------

1.0	06.12.2023	- Creation
2.0	01.02.2024	Update chapters 6, 9, 12, 13

Annex I

Lynk CUFF Short - Article	
02011001	Lynk™ Cuff short, small
02011002	Lynk™ Cuff short, medium
02011003	Lynk™ Cuff short, large

Lynk CUFF Standard - Article	
02012001	Lynk™ Cuff standard, small
02012002	Lynk™ Cuff standard, medium
02012003	Lynk™ Cuff standard, large

Lynk COVER - Article	
02021001	Lynk™ Cover, small
02021002	Lynk™ Cover, medium
02021003	Lynk™ Cover, large

Lynk HOOK Module Standard – Article	
02041000	Lynk™ HOOK Module Standard

Lynk HOLD Module Standard – Article	
02051000	Lynk™ HOLD Module, Standard

Lynk PUSH Module Standard – Article	
02061000	Lynk™ PUSH Module, Standard

Lynk CLIC Module Standard – Article	
02071000	Lynk™ CLIC Module, Standard
02072000	Lynk™ CLIC Head 360°
02073000	Lynk™ CLIC Head 70°
02074000	Lynk™ CLIC Head 0°
02075000	Lynk™ CLIC Adapter Standard

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- Terugsturen naar macu4 GmbH