



# The Explorer System

Order Form



## EXPLORER | SOCKET

Colour

<input type="checkbox"/>  white   yellow	<input type="checkbox"/>  white   red	<input type="checkbox"/>  black   red
<input type="checkbox"/>  black   yellow	<b>Arm side</b> <input type="checkbox"/> right <input type="checkbox"/> left	<b>Thread</b> <input type="checkbox"/> metric <input type="checkbox"/> inch










## EXPLORER | RING

Colour

<input type="checkbox"/>  red	<input type="checkbox"/>  black	<input type="checkbox"/>  yellow
<b>Size</b> <input type="checkbox"/> for Explorer socket <input type="checkbox"/> for third party socket	<b>Thread</b> <input type="checkbox"/> metric <input type="checkbox"/> inch	

## EXPLORER | BIKE MODULE

Colour

<input type="checkbox"/>  black   R-B	<input type="checkbox"/>  black   R-R	<input type="checkbox"/>  black   Y-B
<input type="checkbox"/>  black   Y-Y	<input type="checkbox"/>  red   B-B	<input type="checkbox"/>  red   Y-R
<input type="checkbox"/>  red   Y-Y	<input type="checkbox"/>  blue   R-R	<input type="checkbox"/>  blue   Y-Y
<b>Size</b> <input type="checkbox"/> standard	<b>Thread</b> <input type="checkbox"/> metric <input type="checkbox"/> inch	<b>Spacer</b> ____ (mm)

## EXPLORER | TWIN MODULE

Colour

<input type="checkbox"/>  black   R-B	<input type="checkbox"/>  black   R-R	<input type="checkbox"/>  black   Y-B
<input type="checkbox"/>  black   Y-Y	<input type="checkbox"/>  red   B-B	<input type="checkbox"/>  red   Y-R
<input type="checkbox"/>  red   Y-Y	<input type="checkbox"/>  blue   R-R	<input type="checkbox"/>  blue   Y-Y
<b>Size*</b> <input type="checkbox"/> standard <input type="checkbox"/> large	<b>Thread</b> <input type="checkbox"/> metric <input type="checkbox"/> inch	<b>Spacer</b> ____ (mm)

\* The size should be selected with regard to the handle or bar circumference. See the [Quick Guide for Modules](#).


# The Explorer System

Order Form




## EXPLORER | BALL MODULE


Colour




Size\*\*



Thread



Spacer




## EXPLORER | SWIM MODULE


Colour




Size\*\*



Thread



Spacer



## EXPLORER | PINCH MODULE


Colour




Size



Thread



Spacer\*\*\*



## NOTE

The configuration of the socket and the determination of the spacer length of a module are based on measurement data. Please provide these in the required formats.

If you are using Explorer modules with a third-party socket (not Explorer socket), please enter the required spacer length for each module. The [Quick Guide for third-party sockets](#) provides information on how to determine the spacer length.

Make sure that you order the same thread type for all components so that they are compatible with each other.

Follow the recommendations for using the two Explorer Ring variants. Ring Standard is not compatible with the Explorer socket.

\*\* The size should be selected with regard to the size of the palm. See the [Quick Guide for Modules](#). \*\*\* We recommend not using a spacer if the arm length of the affected arm is less than 50% of the unaffected arm.