

Quick Guide | Module

The Lynk System



Clic Module

- Enables quick attachment and detachment to and from handles / bars
- Available in one size
- With 3 different module Heads and Lamellas



Use it for activities such as:

- Cycling
- Scooter
- Use a broom
- Mini golf, indoor hockey

Hold Module

- Enables smaller and lighter objects to be clamped between the module clamp
- Available in one size
- With integrated ball joint and quick release



Use it for activities such as:

- Use eating cutlery
- Use cooking cutlery
- Do handicrafts
- Rope jumping

Hook Module

- Supports activities where something needs to be carried, pulled or pushed away
- Force transmission to the forearm, not the wrist
- Available in one size



Use it for activities such as:

- Carry a bag
- Pull a suitcase
- Train in the gym

Push Module

- Support for grip fixation and pressure application
- Force transmission to the forearm, not the wrist
- Available in one size
- Fitted with Velcro fastener.



Use it for activities such as:

- Walk with crutches
- Arch shooting, Train in the Gym
- Go kayaking
- Rope jumping